

Annexin IV Rabbit mAb

Catalog No: #59808

Package Size: #59808-1 50ul #59808-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Annexin IV Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	Annexin IV Antibody detects endogenous levels of total Annexin IV
Immunogen Description	A synthesized peptide derived from human Annexin IV
Other Names	AIV; AnnexinA4; AnnexinIV; ANX 4; ANX A4; ANX4; ANXA4; Chromobindin4; EndonexinI; LipocortinIV; P32.5; P33/41; PAPII; PIG28; PP4X; Protein II; Xanx-4; ZAP36;
Accession No.	Uniprot:P09525
Uniprot	P09525
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Annexin IV expression in HepG2 cell lysate.

Product Description

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Background

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Note: This product is for in vitro research use only