

OS9 Rabbit mAb

Catalog No: #59813

Package Size: #59813-1 50ul #59813-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OS9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	OS9 Antibody detects endogenous levels of total OS9
Immunogen Description	A synthesized peptide derived from human OS9
Other Names	ERLEC2; Os9;
Accession No.	Uniprot:Q13438
Uniprot	Q13438
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of OS9 expression in HeLa cell lysate.

Product Description

Lectin which functions in endoplasmic reticulum (ER) quality control and ER-associated degradation (ERAD).

Background

Lectin which functions in endoplasmic reticulum (ER) quality control and ER-associated degradation (ERAD).

Note: This product is for in vitro research use only