

Galectin 7 Rabbit mAb

Catalog No: #59814

Package Size: #59814-1 50ul #59814-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Galectin 7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	Galectin 7 Antibody detects endogenous levels of total Galectin 7
Immunogen Description	A synthesized peptide derived from human Galectin 7
Other Names	Galectin-7; Gal-7; HKL-14; PI7; LGALS7; PIG1;
Accession No.	Uniprot:P47929
Uniprot	P47929
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:100~1:500 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of Galectin 7 expression in HACAT cell lysate.

Product Description

Could be involved in cell-cell and/or cell-matrix interactions necessary for normal growth control. Pro-apoptotic protein that functions intracellularly upstream of JNK activation and cytochrome c release.

Background

Could be involved in cell-cell and/or cell-matrix interactions necessary for normal growth control. Pro-apoptotic protein that functions intracellularly upstream of JNK activation and cytochrome c release.

Note: This product is for in vitro research use only