

CD3G Rabbit mAb

Catalog No: #59821

Package Size: #59821-1 50ul #59821-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD3G Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse
Specificity	CD3G Antibody detects endogenous levels of total CD3G
Immunogen Description	A synthesized peptide derived from human CD3G
Other Names	CD3 gamma; CD3g; T3G;
Accession No.	Uniprot:P09693
Uniprot	P09693
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:100~1:500 ICC/IF 1:100~1:500 IP 1:50 FC 1:60

Images



Western blot analysis of CD3G expression in Jurkat cell lysate.

Product Description

The CD3 complex mediates signal transduction.

Background

The CD3 complex mediates signal transduction.

Note: This product is for in vitro research use only