

NOXA2/p67phox Rabbit mAb

Catalog No: #59828

Package Size: #59828-1 50ul #59828-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NOXA2/p67phox Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	NOXA2/p67phox Antibody detects endogenous levels of total NOXA2/p67phox
Immunogen Description	A synthesized peptide derived from human NOXA2/p67phox
Other Names	Ncf2; NOXA2; P67 PHOX;
Accession No.	Uniprot:P19878
Uniprot	P19878
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:100~1:500 ICC/IF 1:100~1:500 IP 1:50

Images



Western blot analysis of NOXA2/p67phox expression in THP-1 cell lysate.

Product Description

NCF2, NCF1, and a membrane bound cytochrome b558 are required for activation of the latent NADPH oxidase (necessary for superoxide production).

Background

NCF2, NCF1, and a membrane bound cytochrome b558 are required for activation of the latent NADPH oxidase (necessary for superoxide production).

Note: This product is for in vitro research use only