

DAP12 Rabbit mAb

Catalog No: #59832

Package Size: #59832-1 50ul #59832-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DAP12 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	DAP12 Antibody detects endogenous levels of total DAP12
Immunogen Description	A synthesized peptide derived from human DAP12
Other Names	DAP12; KARAP; PLOSL; TYROBP;
Accession No.	Uniprot:O43914
Uniprot	O43914
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:80

Images



Western blot analysis of DAP12 expression in Human PBMC cell lysate.

Product Description

Non-covalently associates with activating receptors of the CD300 family. Cross-linking of CD300-TYROBP complexes results in cellular activation.

Background

Non-covalently associates with activating receptors of the CD300 family. Cross-linking of CD300-TYROBP complexes results in cellular activation.

Note: This product is for in vitro research use only