

## ORP150 Rabbit mAb

Catalog No: #59844

Package Size: #59844-1 50ul #59844-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ORP150 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Mouse Rat
Specificity	ORP150 Antibody detects endogenous levels of total ORP150
Immunogen Description	A synthesized peptide derived from human ORP150
Other Names	Grp170; HSP12A; Hyou1; Orp150;
Accession No.	Uniprot:Q9Y4L1
Uniprot	Q9Y4L1
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 FC 1:50

## Images



Western blot analysis of ORP150 expression in MCF7 cell lysate.

## Product Description

Has a pivotal role in cytoprotective cellular mechanisms triggered by oxygen deprivation. May play a role as a molecular chaperone and participate in protein folding.

## Background

---

Has a pivotal role in cytoprotective cellular mechanisms triggered by oxygen deprivation. May play a role as a molecular chaperone and participate in protein folding.

---

Note: This product is for in vitro research use only