

ADAMTS13 Rabbit mAb

Catalog No: #59847

Package Size: #59847-1 50ul #59847-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ADAMTS13 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	ADAMTS13 Antibody detects endogenous levels of total ADAMTS13
Immunogen Description	A synthesized peptide derived from human ADAMTS13
Other Names	ADAMTS13; ADAM-TS; ADAM-TS13; ADAMTS-13; ADAM-TS 13; C9orf8; TTP; VWF-CP; VWFCP;
Accession No.	Uniprot:Q76LX8
Uniprot	Q76LX8
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:100

Images



Western blot analysis of ADAMTS13 expression in A549 cell lysate.

Product Description

Cleaves the vWF multimers in plasma into smaller forms.

Background

Cleaves the vWF multimers in plasma into smaller forms.

Note: This product is for in vitro research use only