

CD35 Rabbit mAb

Catalog No: #59855

Package Size: #59855-1 50ul #59855-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD35 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	CD35 Antibody detects endogenous levels of total CD35
Immunogen Description	A synthesized peptide derived from human CD35
Other Names	C3 binding protein; C3b/C4b receptor; C3BR; C4BR; CD35; CR1; KN;
Accession No.	Uniprot:P17927
Uniprot	P17927
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of CD35 expression in TF1 cell lysate.

Product Description

Mediates cellular binding of particles and immune complexes that have activated complement.

Background

Mediates cellular binding of particles and immune complexes that have activated complement.

Note: This product is for in vitro research use only