

FKBP51 Rabbit mAb

Catalog No: #59856

Package Size: #59856-1 50ul #59856-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FKBP51 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Rat
Specificity	FKBP51 Antibody detects endogenous levels of total FKBP51
Immunogen Description	A synthesized peptide derived from human FKBP51
Other Names	AIG6; FKBP5; FKBP51; FKBP54; P54; PPlase; PPlase FKBP5; Ptg10;Rotamase;
Accession No.	Uniprot:Q13451
Uniprot	Q13451
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of FKBP51 expression in Caco-2 cell lysate.

Product Description

Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.

Background

Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.

Note: This product is for in vitro research use only