

FKBP52 Rabbit mAb

Catalog No: #59857

Package Size: #59857-1 50ul #59857-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FKBP52 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	FKBP52 Antibody detects endogenous levels of total FKBP52
Immunogen Description	A synthesized peptide derived from human FKBP52
Other Names	HBI; p52; Hsp56; FKBP51; FKBP52; FKBP59; PPIase;
Accession No.	Uniprot:Q02790
Uniprot	Q02790
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of FKBP52 expression in MCF7 cell lysate.

Product Description

Immunophilin protein with PPIase and co-chaperone activities (By similarity). Component of unliganded steroid receptors heterocomplexes through interaction with heat-shock protein 90 (HSP90).

Background

Immunophilin protein with PPIase and co-chaperone activities (By similarity). Component of unliganded steroid receptors heterocomplexes through interaction with heat-shock protein 90 (HSP90).

Note: This product is for in vitro research use only