

ADAR1 Rabbit mAb

Catalog No: #59865

Package Size: #59865-1 50ul #59865-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ADAR1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	ADAR1 Antibody detects endogenous levels of total ADAR1
Immunogen Description	A synthesized peptide derived from human ADAR1
Other Names	ADAR; Adar1; AGS6; DRADA; Dsh; Dsrad; IFI4; P136;
Accession No.	Uniprot:P55265
Uniprot	P55265
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of ADAR1 expression in Ramos cell lysate.

Product Description

Converts multiple adenosines to inosines and creates I/U mismatched base pairs in double-helical RNA substrates without apparent sequence specificity.

Background

Converts multiple adenosines to inosines and creates I/U mismatched base pairs in double-helical RNA substrates without apparent sequence specificity.

Note: This product is for in vitro research use only