

MGEA5 Rabbit mAb

Catalog No: #59869

Package Size: #59869-1 50ul #59869-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MGEA5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse
Specificity	MGEA5 Antibody detects endogenous levels of total MGEA5
Immunogen Description	A synthesized peptide derived from human MGEA5
Other Names	HEXC; MEA5; Mgea5; NCOAT; OGA;
Accession No.	Uniprot:O60502
Uniprot	O60502
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 IP 1:80

Images



Western blot analysis of MGEA5 expression in JAR cell lysate.

Product Description

Isoform 1: Cleaves GlcNAc but not GalNAc from O-glycosylated proteins. Isoform 3: Cleaves GlcNAc but not GalNAc from O-glycosylated proteins.

Background

Isoform 1: Cleaves GlcNAc but not GalNAc from O-glycosylated proteins. Isoform 3: Cleaves GlcNAc but not GalNAc from O-glycosylated proteins.

Note: This product is for in vitro research use only