

## Frizzled 8 Rabbit mAb

Catalog No: #59872

Package Size: #59872-1 50ul #59872-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Frizzled 8 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	Frizzled 8 Antibody detects endogenous levels of total Frizzled 8
Immunogen Description	A synthesized peptide derived from human Frizzled 8
Other Names	Frizzled-8; FZ8; FZD8; hFz8;
Accession No.	Uniprot:Q9H461
Uniprot	Q9H461
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200

## Images



Western blot analysis of Frizzled 8 expression in Jurkat cell lysate.

## Product Description

Receptor for Wnt proteins. Component of the Wnt-Fzd-LRP5-LRP6 complex that triggers beta-catenin signaling through inducing aggregation of receptor-ligand complexes into ribosome-sized signalosomes.

## Background

---

Receptor for Wnt proteins. Component of the Wnt-Fzd-LRP5-LRP6 complex that triggers beta-catenin signaling through inducing aggregation of receptor-ligand complexes into ribosome-sized signalosomes.

---

Note: This product is for in vitro research use only