

BRD2 Rabbit mAb

Catalog No: #59880

Package Size: #59880-1 50ul #59880-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	BRD2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	BRD2 Antibody detects endogenous levels of total BRD2
Immunogen Description	A synthesized peptide derived from human BRD2
Other Names	Brd2; FSH; FSRG1; NAT; O27.1.; RING3; RNF3;
Accession No.	Uniprot:P25440
Uniprot	P25440
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:100

Images



Western blot analysis of BRD2 expression in NCCIT cell lysate.

Product Description

May play a role in spermatogenesis or folliculogenesis.

Background

May play a role in spermatogenesis or folliculogenesis.

Note: This product is for in vitro research use only