

Coilin Rabbit mAb

Catalog No: #59884

Package Size: #59884-1 50ul #59884-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Coilin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	Coilin Antibody detects endogenous levels of total Coilin
Immunogen Description	A synthesized peptide derived from human Coilin
Other Names	CLN80; COIL; Coilin; Coilin p80;
Accession No.	Uniprot:P38432
Uniprot	P38432
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:100

Images



Western blot analysis of Coilin expression in Jurkat cell lysate.

Product Description

Is a component of the nuclear coiled bodies (CBS) which are involved in the function or assembly/disassembly of nucleoplasmic snRNPs. During mitosis, CBS disassemble, coinciding with a mitotic-specific phosphorylation of p80 coilin.

Background

Is a component of the nuclear coiled bodies (CBS) which are involved in the function or assembly/disassembly of nucleoplasmic snRNPs. During mitosis, CBS disassemble, coinciding with a mitotic-specific phosphorylation of p80 coilin.

Note: This product is for in vitro research use only