

## PGAM1 Rabbit mAb

Catalog No: #59889

Package Size: #59889-1 50ul #59889-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PGAM1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	PGAM1 Antibody detects endogenous levels of total PGAM1
Immunogen Description	A synthesized peptide derived from human PGAM1
Other Names	PGAM1; PGAMA; PGAMB;
Accession No.	Uniprot:P18669
Uniprot	P18669
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 FC 1:100

## Images



Western blot analysis of PGAM1 expression in A431 cell lysate.

## Product Description

Interconversion of 3- and 2-phosphoglycerate with 2,3-bisphosphoglycerate as the primer of the reaction. Can also catalyze the reaction of EC 5.4.2.4 (synthase) and EC 3.1.3.13 (phosphatase), but with a reduced activity.

## Background

---

Interconversion of 3- and 2-phosphoglycerate with 2,3-bisphosphoglycerate as the primer of the reaction. Can also catalyze the reaction of EC 5.4.2.4 (synthase) and EC 3.1.3.13 (phosphatase), but with a reduced activity.

---

Note: This product is for in vitro research use only