

PRKAR1A Rabbit mAb

Catalog No: #59898

Package Size: #59898-1 50ul #59898-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PRKAR1A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	PRKAR1A Antibody detects endogenous levels of total PRKAR1A
Immunogen Description	A synthesized peptide derived from human PRKAR1A
Other Names	PRKAR1A; CAR; CNC; CNC1; PKR1; PPNAD1; PRKAR1; TSE1;
Accession No.	Uniprot:P10644
Uniprot	P10644
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of Protein Kinase A regulatory subunit I alpha expression in U-87MG cell lysate.

Product Description

Regulatory subunit of the cAMP-dependent protein kinases involved in cAMP signaling in cells.

Background

Regulatory subunit of the cAMP-dependent protein kinases involved in cAMP signaling in cells.

Note: This product is for in vitro research use only