

Syntaxin 3 Rabbit mAb

Catalog No: #59899

Package Size: #59899-1 50ul #59899-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Syntaxin 3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Syntaxin 3 Antibody detects endogenous levels of total Syntaxin 3
Immunogen Description	A synthesized peptide derived from human Syntaxin 3
Other Names	Stx3; STX3A; Syntaxin 3A;
Accession No.	Uniprot:Q13277
Uniprot	Q13277
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Syntaxin 3 expression in HepG2 cell lysate.

Product Description

Potentially involved in docking of synaptic vesicles at presynaptic active zones.

Background

Potentially involved in docking of synaptic vesicles at presynaptic active zones.

Note: This product is for in vitro research use only