

CD27 Rabbit mAb

Catalog No: #59901

Package Size: #59901-1 50ul #59901-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD27 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Rat
Specificity	CD27 Antibody detects endogenous levels of total CD27
Immunogen Description	A synthesized peptide derived from human CD27
Other Names	CD27; CD27L receptor; LPFS2; S152; T14; TNFRSF7; TNFSF7; Tp55;
Accession No.	Uniprot:P26842
Uniprot	P26842
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of CD27 expression in Ramos cell lysate.

Product Description

Receptor for CD70/CD27L. May play a role in survival of activated T-cells. May play a role in apoptosis through association with SIVA1.

Background

Receptor for CD70/CD27L. May play a role in survival of activated T-cells. May play a role in apoptosis through association with SIVA1.

Note: This product is for in vitro research use only