

SFRP4 Rabbit mAb

Catalog No: #59911

Package Size: #59911-1 50ul #59911-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SFRP4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	SFRP4 Antibody detects endogenous levels of total SFRP4
Immunogen Description	A synthesized peptide derived from human SFRP4
Other Names	FRP4; FrpHE; PYL; Sfrp4;
Accession No.	Uniprot:Q6FHJ7
Uniprot	Q6FHJ7
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of SFRP4 expression in 293T cell lysate.

Product Description

Soluble frizzled-related proteins (sFRPS) function as modulators of Wnt signaling through direct interaction with Wnts. They have a role in regulating cell growth and differentiation in specific cell types.

Background

Soluble frizzled-related proteins (sFRPS) function as modulators of Wnt signaling through direct interaction with Wnts. They have a role in regulating cell growth and differentiation in specific cell types.

Note: This product is for in vitro research use only