

GAB1 Rabbit mAb

Catalog No: #59919

Package Size: #59919-1 50ul #59919-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GAB1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	GAB1 Antibody detects endogenous levels of total GAB1
Immunogen Description	A synthesized peptide derived from human GAB1
Other Names	GAB1; GRB 2 associated binder 1;
Accession No.	Uniprot:Q13480
Uniprot	Q13480
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:500

Images



Western blot analysis of GAB1 expression in HEK293 cell lysate.

Product Description

Probably involved in EGF and insulin receptor signaling.

Background

Probably involved in EGF and insulin receptor signaling.

Note: This product is for in vitro research use only