

CA125 / MUC16 Rabbit mAb

Catalog No: #59921

Package Size: #59921-1 50ul #59921-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CA125 / MUC16 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	CA-125 (MUC16) Antibody detects endogenous levels of total CA-125 (MUC16)
Immunogen Description	A synthesized peptide derived from human CA-125 (MUC16)
Other Names	CA125; MUC16; Mucin 16; Mucin16;
Accession No.	Uniprot:Q8WXI7
Uniprot	Q8WXI7
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Immunohistochemical analysis of paraffin-embedded human ovary cancer, using CA-125 (MUC16) Antibody.

Product Description

Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces.

Background

Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces.

Note: This product is for in vitro research use only