

FOXP1 Rabbit mAb

Catalog No: #59922

Package Size: #59922-1 50ul #59922-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FOXP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	FOXP1 Antibody detects endogenous levels of total FOXP1
Immunogen Description	A synthesized peptide derived from human FOXP1
Other Names	FOXP 1; Glutamine rich factor 1; hFKH1B; HSPC215; QRF1;
Accession No.	Uniprot:Q9H334
Uniprot	Q9H334
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:100~1:500 IP 1:50 FC:1:50

Images



Western blot analysis of FOXP1 expression in Daudi cell lysate.

Product Description

Transcriptional repressor. It plays an important role in the specification and differentiation of lung epithelium. Can act with CTBP1 to synergistically repress transcription but CTPBP1 is not essential. Essential transcriptional regulator of B cell development.

Background

Transcriptional repressor. It plays an important role in the specification and differentiation of lung epithelium. Can act with CTBP1 to synergistically repress transcription but CTPBP1 is not essential. Essential transcriptional regulator of B cell development.

Note: This product is for in vitro research use only