

Villin-1 Rabbit mAb

Catalog No: #59924

Package Size: #59924-1 50ul #59924-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Villin-1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	Villin-1 Antibody detects endogenous levels of total Villin-1
Immunogen Description	A synthesized peptide derived from human Villin-1
Other Names	VIL; VIL1; Villin1;
Accession No.	Uniprot:P09327
Uniprot	P09327
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:10000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of Villin-1 expression in CaCo-2 cell lysate.

Product Description

Ca(2+)-regulated actin-binding protein.

Background

Ca(2+)-regulated actin-binding protein.

Note: This product is for in vitro research use only