

Alpha actinin 2 Rabbit mAb

Catalog No: #59941

Package Size: #59941-1 50ul #59941-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Alpha actinin 2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Mouse, Rat, Human
Specificity	Alpha actinin 2 Antibody detects endogenous levels of total Alpha actinin 2
Immunogen Description	A synthesized peptide derived from human Alpha actinin 2
Other Names	Actin binding protein; Actinin alpha 2; ACTN2; Alpha actinin 2; CMD1AA;
Accession No.	Uniprot:P35609
Uniprot	P35609
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 IP 1:30

Images



Western blot analysis of Alpha actinin 2 expression in human skeletal muscle lysate.

Product Description

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Background

F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

Note: This product is for in vitro research use only