

HB EGF Rabbit mAb

Catalog No: #59942

Package Size: #59942-1 50ul #59942-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HB EGF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Rat
Specificity	HB EGF Antibody detects endogenous levels of total HB EGF
Immunogen Description	A synthesized peptide derived from human HB EGF
Other Names	DTR; DTS; DTSF; HB-EGF; HBEGF; HEGFL;
Accession No.	Uniprot:Q99075
Uniprot	Q99075
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of HB EGF expression in K562 cell lysate.

Product Description

May be involved in macrophage-mediated cellular proliferation. It is mitogenic for fibroblasts and smooth muscle but not endothelial cells.

Background

May be involved in macrophage-mediated cellular proliferation. It is mitogenic for fibroblasts and smooth muscle but not endothelial cells.

Note: This product is for in vitro research use only