

USP22 Rabbit mAb

Catalog No: #59947

Package Size: #59947-1 50ul #59947-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	USP22 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	USP22 Antibody detects endogenous levels of total USP22
Immunogen Description	A synthesized peptide derived from human USP22
Other Names	Ubiquitin carboxyl terminal hydrolase 22; Usp22; USP3L;
Accession No.	Uniprot:Q9UPT9
Uniprot	Q9UPT9
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of USP22 expression in (1) HeLa cell lysate; (2) Mouse spleen lysate.

Product Description

Histone deubiquitinating component of the transcription regulatory histone acetylation (HAT) complex SAGA. Catalyzes the deubiquitination of both histones H2A and H2B, thereby acting as a coactivator.

Background

Histone deubiquitinating component of the transcription regulatory histone acetylation (HAT) complex SAGA. Catalyzes the deubiquitination of both histones H2A and H2B, thereby acting as a coactivator.

Note: This product is for in vitro research use only