

Jarid2 Rabbit mAb

Catalog No: #59952

Package Size: #59952-1 50ul #59952-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Jarid2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	Jarid2 Antibody detects endogenous levels of total Jarid2
Immunogen Description	A synthesized peptide derived from human Jarid2
Other Names	JARD2; JARID2; JMJ;
Accession No.	Uniprot:Q92833
Uniprot	Q92833
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

Images



Western blot analysis of Jarid2 expression in NCCIT cell lysate.

Product Description

Regulator of histone methyltransferase complexes that plays an essential role in embryonic development, including heart and liver development, neural tube fusion process and hematopoiesis.

Background

Regulator of histone methyltransferase complexes that plays an essential role in embryonic development, including heart and liver development, neural tube fusion process and hematopoiesis.

Note: This product is for in vitro research use only