

## Monoamine Oxidase A Rabbit mAb

Catalog No: #59953

Package Size: #59953-1 50ul #59953-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Monoamine Oxidase A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Monoamine Oxidase A Antibody detects endogenous levels of total Monoamine Oxidase A
Immunogen Description	A synthesized peptide derived from human Monoamine Oxidase A
Other Names	AOFA; maoA;
Accession No.	Uniprot:P21397
Uniprot	P21397
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

## Images



Western blot analysis of Monoamine Oxidase A expression in HepG2 cell lysate.

## Product Description

Catalyzes the oxidative deamination of biogenic and xenobiotic amines and has important functions in the metabolism of neuroactive and vasoactive amines in the central nervous system and peripheral tissues.

## Background

---

Catalyzes the oxidative deamination of biogenic and xenobiotic amines and has important functions in the metabolism of neuroactive and vasoactive amines in the central nervous system and peripheral tissues.

---

Note: This product is for in vitro research use only