

IFNGR1 Rabbit mAb

Catalog No: #59954

Package Size: #59954-1 50ul #59954-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IFNGR1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	IFNGR1 Antibody detects endogenous levels of total IFNGR1
Immunogen Description	A synthesized peptide derived from human IFNGR1
Other Names	CD119; CDw119; IFN gamma R alpha; IFNGR1;
Accession No.	Uniprot:P15260
Uniprot	P15260
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:500

Images



Western blot analysis of IFNGR1 expression in MCF7 cell lysate.

Product Description

Receptor for interferon gamma. Two receptors bind one interferon gamma dimer.

Background

Receptor for interferon gamma. Two receptors bind one interferon gamma dimer.

Note: This product is for in vitro research use only