

ICAD Rabbit mAb

Catalog No: #59967

Package Size: #59967-1 50ul #59967-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ICAD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	ICAD Antibody detects endogenous levels of total ICAD
Immunogen Description	A synthesized peptide derived from human ICAD
Other Names	DFF1; DFF35; DFF45; DFFA; H13; ICAD; ICAD L; ICAD S;
Accession No.	Uniprot:O00273
Uniprot	O00273
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of ICAD expression in HeLa cell lysate.

Product Description

Inhibitor of the caspase-activated DNase (DFF40).

Background

Inhibitor of the caspase-activated DNase (DFF40).

Note: This product is for in vitro research use only