

Oct-2 Rabbit mAb

Catalog No: #59971

Package Size: #59971-1 50ul #59971-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Oct-2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC FC
Species Reactivity	Human
Specificity	Oct-2 Antibody detects endogenous levels of total Oct-2
Immunogen Description	A synthesized peptide derived from human Oct-2
Other Names	Oct 2; Octamer binding transcription factor 2; OTF2; POU domain; POU2F2;
Accession No.	Uniprot:P09086
Uniprot	P09086
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Oct-2 expression in Daudi cell lysate.

Product Description

Transcription factor that specifically binds to the octamer motif (5'-ATTTGCAT-3'). Regulates transcription in a number of tissues in addition to activating immunoglobulin gene expression. Modulates transcription transactivation by NR3C1, AR and PGR.

Background

Transcription factor that specifically binds to the octamer motif (5'-ATTGTCAT-3'). Regulates transcription in a number of tissues in addition to activating immunoglobulin gene expression. Modulates transcription transactivation by NR3C1, AR and PGR.

Note: This product is for in vitro research use only