

Activin Receptor Type IIB/ACVR2B Rabbit mAb

Catalog No: #59980

Package Size: #59980-1 50ul #59980-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Activin Receptor Type IIB/ACVR2B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Activin Receptor Type IIB Antibody detects endogenous levels of total Activin Receptor Type IIB
Immunogen Description	A synthesized peptide derived from human Activin Receptor Type IIB
Other Names	ActR IIB; ACVR 2B; HTX4;
Accession No.	Uniprot:Q13705
Uniprot	Q13705
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Activin Receptor Type IIB/ACVR2B expression in (1) K562 cell lysate; (2) Mouse heart lysate.

Product Description

On ligand binding, forms a receptor complex consisting of two type II and two type I transmembrane serine/threonine kinases. Type II receptors phosphorylate and activate type I receptors which autophosphorylate, then bind and activate SMAD transcriptional regulators. Receptor for activin A, activin B and inhibin A.

Background

On ligand binding, forms a receptor complex consisting of two type II and two type I transmembrane serine/threonine kinases. Type II receptors phosphorylate and activate type I receptors which autophosphorylate, then bind and activate SMAD transcriptional regulators. Receptor for activin A, activin B and inhibin A.

Note: This product is for in vitro research use only