

## IGFBP1 Rabbit mAb

Catalog No: #59986

Package Size: #59986-1 50ul #59986-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	IGFBP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	IGFBP1 Antibody detects endogenous levels of total IGFBP1
Immunogen Description	A synthesized peptide derived from human IGFBP1
Other Names	AFBP; IBP1; PP12; IGF-BP25; hIGFBP-1;
Accession No.	Uniprot:P08833
Uniprot	P08833
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of IGFBP1 expression in HepG2 cell lysate.

## Product Description

IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors. Promotes cell migration.

## Background

---

IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors. Promotes cell migration.

---

Note: This product is for in vitro research use only