

IL7R alpha/CD127 Rabbit mAb

Catalog No: #59988

Package Size: #59988-1 50ul #59988-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | IL7R alpha/CD127 Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB FC |
| Species Reactivity | Human Mouse Rat |
| Specificity | IL7R alpha/CD127 Antibody detects endogenous levels of total IL7R alpha/CD127 |
| Immunogen Description | A synthesized peptide derived from human IL7R alpha/CD127 |
| Other Names | IL7R; CD127; CDW127; IL-7R-alpha; IL7RA; ILRA ; |
| Accession No. | Uniprot:P16871 |
| Uniprot | P16871 |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

WB 1:500~1:2000 FC 1:100

Images



Western blot analysis of IL7R alpha/CD127 expression in (1) HeLa cell lysate; (2) Mouse spleen lysate.

Product Description

Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

Background

Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

Note: This product is for in vitro research use only