

## IL7R alpha/CD127 Rabbit mAb

Catalog No: #59988

Package Size: #59988-1 50ul #59988-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	IL7R alpha/CD127 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	IL7R alpha/CD127 Antibody detects endogenous levels of total IL7R alpha/CD127
Immunogen Description	A synthesized peptide derived from human IL7R alpha/CD127
Other Names	IL7R; CD127; CDW127; IL-7R-alpha; IL7RA; ILRA ;
Accession No.	Uniprot:P16871
Uniprot	P16871
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 FC 1:100

## Images



Western blot analysis of IL7R alpha/CD127 expression in (1) HeLa cell lysate; (2) Mouse spleen lysate.

## Product Description

Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

## Background

Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

---

Note: This product is for in vitro research use only