

KAISO Rabbit mAb

Catalog No: #60003

Package Size: #60003-1 50ul #60003-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	KAISO Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	KAISO Antibody detects endogenous levels of total KAISO
Immunogen Description	A synthesized peptide derived from human KAISO
Other Names	Kaiso; Zbtb33; ZNF Kaiso; ZNF348;
Accession No.	Uniprot:Q86T24
Uniprot	Q86T24
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of KAISO expression in MCF7 cell lysate.

Product Description

Transcriptional regulator with bimodal DNA-binding specificity. Binds to methylated CpG dinucleotides in the consensus sequence 5'-CGCG-3' and also binds to the non-methylated consensus sequence 5'-CTGCNA-3'.

Background

Transcriptional regulator with bimodal DNA-binding specificity. Binds to methylated CpG dinucleotides in the consensus sequence 5'-CGCG-3' and also binds to the non-methylated consensus sequence 5'-CTGCNA-3'.

Note: This product is for in vitro research use only