

CRKL Rabbit mAb

Catalog No: #60007

Package Size: #60007-1 50ul #60007-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CRKL Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	CRKL Antibody detects endogenous levels of total CRKL
Immunogen Description	A synthesized peptide derived from human CRKL
Other Names	Crk L; Crk like protein; Crkl; Crkol;
Accession No.	Uniprot:P46109
Uniprot	P46109
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of CRKL expression in K562 cell lysate.

Product Description

May mediate the transduction of intracellular signals.

Background

May mediate the transduction of intracellular signals.

Note: This product is for in vitro research use only