

## HRH3 Rabbit mAb

Catalog No: #60014

Package Size: #60014-1 50ul #60014-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HRH3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	HRH3 Antibody detects endogenous levels of total HRH3
Immunogen Description	A synthesized peptide derived from human HRH3
Other Names	GPCR97; H3R; HH3R; Histamine H3 receptor; HRH3;
Accession No.	Uniprot:Q9Y5N1
Uniprot	Q9Y5N1
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:100

## Images



Western blot analysis of HRH3 expression in A549 cell lysate.

## Product Description

The H3 subclass of histamine receptors could mediate the histamine signals in CNS and peripheral nervous system. Signals through the inhibition of adenylate cyclase and displays high constitutive activity (spontaneous activity in the absence of agonist).

## Background

---

The H3 subclass of histamine receptors could mediate the histamine signals in CNS and peripheral nervous system. Signals through the inhibition of adenylate cyclase and displays high constitutive activity (spontaneous activity in the absence of agonist).

---

Note: This product is for in vitro research use only