

ZNF217 Rabbit mAb

Catalog No: #60017

Package Size: #60017-1 50ul #60017-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | ZNF217 Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB |
| Species Reactivity | Human |
| Specificity | ZNF217 Antibody detects endogenous levels of total ZNF217 |
| Immunogen Description | A synthesized peptide derived from human ZNF217 |
| Other Names | ZABC1; ZN217; ZNF217; |
| Accession No. | Uniprot:O75362 |
| Uniprot | O75362 |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

WB 1:500~1:2000

Images



Western blot analysis of ZNF217 expression in Caco2 cell lysate.

Product Description

Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death. Promotes phosphorylation of AKT1 at 'Ser-473'.

Background

Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death. Promotes phosphorylation of AKT1 at 'Ser-473'.

Note: This product is for in vitro research use only