

ZNF217 Rabbit mAb

Catalog No: #60017

Package Size: #60017-1 50ul #60017-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF217 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	ZNF217 Antibody detects endogenous levels of total ZNF217
Immunogen Description	A synthesized peptide derived from human ZNF217
Other Names	ZABC1; ZN217; ZNF217;
Accession No.	Uniprot:O75362
Uniprot	O75362
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of ZNF217 expression in Caco2 cell lysate.

Product Description

Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death. Promotes phosphorylation of AKT1 at 'Ser-473'.

Background

Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death. Promotes phosphorylation of AKT1 at 'Ser-473'.

Note: This product is for in vitro research use only