

Olig1 Rabbit mAb

Catalog No: #60018

Package Size: #60018-1 50ul #60018-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Olig1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	Olig1 Antibody detects endogenous levels of total Olig1
Immunogen Description	A synthesized peptide derived from human Olig1
Other Names	BHLH B6; BHLHB 6; bHLHb6; bHLHe21; Olig1; Oligo1;
Accession No.	Uniprot:Q8TAK6
Uniprot	Q8TAK6
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:50

Images



Western blot analysis of Olig1 expression in MOLT4 cell lysate.

Product Description

Promotes formation and maturation of oligodendrocytes, especially within the brain. Cooperates with OLIG2 to establish the pMN domain of the embryonic neural tube.

Background

Promotes formation and maturation of oligodendrocytes, especially within the brain. Cooperates with OLIG2 to establish the pMN domain of the embryonic neural tube.

Note: This product is for in vitro research use only