

GLB1 Rabbit mAb

Catalog No: #60028

Package Size: #60028-1 50ul #60028-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GLB1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	GLB1 Antibody detects endogenous levels of total GLB1
Immunogen Description	A synthesized peptide derived from human GLB1
Other Names	EBP; ELNR1; Galactosidase beta 1; GLB1; Lactase; MPS4B;
Accession No.	Uniprot:P16278
Uniprot	P16278
Calculated MW	95kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB 1:1000~1:5000 FC 1:50

Images



Western blot analysis of GLB1 expression in SH-SY5Y cell lysate.

All lanes use the Antibody at 1:3K dilution for 1 hour at room temperature.



Product Description

Cleaves beta-linked terminal galactosyl residues from gangliosides, glycoproteins, and glycosaminoglycans.

Background

Cleaves beta-linked terminal galactosyl residues from gangliosides, glycoproteins, and glycosaminoglycans.

Note: This product is for in vitro research use only