

## STK39 Rabbit mAb

Catalog No: #60030

Package Size: #60030-1 50ul #60030-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	STK39 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human
Specificity	STK39 Antibody detects endogenous levels of total STK39
Immunogen Description	A synthesized peptide derived from human STK39
Other Names	DCHT; PASK; SPAK; Stk39;
Accession No.	Uniprot:Q9UEW8
Uniprot	Q9UEW8
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 FC 1:50

## Images



Western blot analysis of STK39 expression in HepG2 cell lysate.

## Product Description

May act as a mediator of stress-activated signals.

## Background

May act as a mediator of stress-activated signals.

---

Note: This product is for in vitro research use only