

PLAP Rabbit mAb

Catalog No: #60031

Package Size: #60031-1 50ul #60031-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PLAP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	PLAP Antibody detects endogenous levels of total PLAP
Immunogen Description	A synthesized peptide derived from human PLAP
Other Names	Glycerophosphatase; PLA2P; Plaa; PLAP;
Accession No.	Uniprot:Q9Y263
Uniprot	Q9Y263
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of PLAP expression in MCF7 cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Involved in the maintenance of ubiquitin levels.

Background

Involved in the maintenance of ubiquitin levels.

Note: This product is for in vitro research use only