

## eIF3B Rabbit mAb

Catalog No: #60032

Package Size: #60032-1 50ul #60032-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	eIF3B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	eIF3B Antibody detects endogenous levels of total eIF3B
Immunogen Description	A synthesized peptide derived from human eIF3B
Other Names	eIF3b; EIF3S9; hPrt1; PRT1;
Accession No.	Uniprot:P55884
Uniprot	P55884
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

## Images



Western blot analysis of eIF3B expression in (1) A431 cell lysate; (2) 293T cell lysate.

## Product Description

Component of the eukaryotic translation initiation factor 3 (eIF-3) complex, which is required for several steps in the initiation of protein synthesis.

## Background

Component of the eukaryotic translation initiation factor 3 (eIF-3) complex, which is required for several steps in the initiation of protein synthesis.

---

Note: This product is for in vitro research use only