

SLC27A4 / FATP4 Rabbit mAb

Catalog No: #60043

Package Size: #60043-1 50ul #60043-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SLC27A4 / FATP4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	SLC27A4 / FATP4 Antibody detects endogenous levels of total SLC27A4 / FATP4
Immunogen Description	A synthesized peptide derived from human SLC27A4 / FATP4
Other Names	ACSVL4; FATP4; Fatty acid transport protein 4; IPS; S27A4; Solute carrier family 27 member4;
Accession No.	Uniprot:Q6P1M0
Uniprot	Q6P1M0
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:100

Images



Western blot analysis of SLC27A4 / FATP4 expression in (1) HeLa cell lysate; (2) Mouse kidney lysate.

Product Description

plays a role in the transport of long chain fatty acids across the plasma membrane. It has acyl-coA ligase activity for long chain and very long chain fatty acids.

Background

plays a role in the transport of long chain fatty acids across the plasma membrane. It has acyl-coA ligase activity for long chain and very long chain fatty acids.

Note: This product is for in vitro research use only