

## TRAF2 Rabbit mAb

Catalog No: #60047

Package Size: #60047-1 50ul #60047-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TRAF2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	TRAF2 Antibody detects endogenous levels of total TRAF2
Immunogen Description	A synthesized peptide derived from human TRAF2
Other Names	TNF receptor associated factor 2; TRAF2; TRAP; TRAP3;
Accession No.	Uniprot:Q12933
Uniprot	Q12933
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

## Images



Western blot analysis of TRAF2 expression in Raji cell lysate.

## Product Description

Regulates activation of NF-kappa-B and JNK and plays a central role in the regulation of cell survival and apoptosis.

## Background

Regulates activation of NF-kappa-B and JNK and plays a central role in the regulation of cell survival and apoptosis.

---

Note: This product is for in vitro research use only