

RCHY1 Rabbit mAb

Catalog No: #60049

Package Size: #60049-1 50ul #60049-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RCHY1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human Mouse
Specificity	RCHY1 Antibody detects endogenous levels of total RCHY1
Immunogen Description	A synthesized peptide derived from human RCHY1
Other Names	ARNIP; CHIMP; hARNIP; hPirh2; PIRH2E; PIRH2F; RCHY1; RNF199; ZCHY; ZNF363;
Accession No.	Uniprot:Q96PM5
Uniprot	Q96PM5
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of RCHY1 expression in 293 cell lysate.

Product Description

Mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including p53/TP53, HDAC1 and CDKN1B. Preferentially acts on tetrameric p53/TP53.

Background

Mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including p53/TP53, HDAC1 and CDKN1B. Preferentially acts on tetrameric p53/TP53.

Note: This product is for in vitro research use only